


Miniatur-Panelmontage-Thermoelement-Steckverbinder mit Edelstahlhalterung JM-E-SSPF Typ E JIS



Miniature in-line Socket with stainless steel mounting bracket
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Panelmontage-Thermoelement-Steckverbinder mit Edelstahlhalterung JM-E-SSPF Typ E JIS

Edelstahlhalterung, die einfach am Verbinder befestigt ist. Vollständig geschlossene Anschlüsse mit Kabelanschlussdichtung

Technische Daten**Specifications**

| | |
|----------------------------|---|
| Product Code | XE-9238-001 |
| General Description | Fully enclosed terminals with cable port seal |
| Thermocouple Type | E JIS |
| Connector Type | Socket with SS Bracket |
| Connector Size | Miniature |
| Colour | Violet |
| Max. Temperature | 220°C |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |