


Standard-Thermoelementstecker JS-E-M Typ E JIS



Standard in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelementstecker JS-E-M Typ E JIS

Die Kontakte bestehen aus echten Thermoelementlegierungen. Geeignet für industrielle Anwendungen

Technische Daten**Specifications**

| | |
|----------------------------|--|
| Product Code | XE-9350-001 |
| General Description | Contacts are polarised to prevent incorrect connection |
| Thermocouple Type | E JIS |
| Connector Type | Plug |
| Connector Size | Standard |
| Colour | Violet |
| Max. Temperature | 220°C |
| Pins | Solid Round Pins |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Copper Nickel |
| Standards Met | BS EN 50212:1997. RoHS |