



Standard-Thermoelement-Steckverbinder In-Line-Buchse JS-J-F Typ J JIS



Standard in-line Socket
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Thermoelement-Steckverbinder In-Line-Buchse JS-J-F Typ J JIS

Akzeptiert standardmige runde Steckerstifte, geeignet fr industrielle Anwendungen

Technische Daten**Specifications**

Product Code	XE-9292-001
General Description	Accepts standard round plug pins, contacts are polarised to prevent incorrect connection.
Thermocouple Type	J JIS
Connector Type	Socket
Connector Size	Standard
Colour	Yellow
Max. Temperature	220°C
+Positive Contact	Iron
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS