


Standard-Rechteckige Faszien-Thermoelement-Steckdose JS-E-FF Typ E JIS



Standard Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Standard-Rechteckiger Thermoelement-Steckverbinder-Faszien-Sockel JS-E-FF Typ E JIS

Nylonclip für schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich große Pins polarisiert

Technische Daten**Specifications**

Product Code	XE-9357-001
General Description	Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel
Thermocouple Type	E JIS
Connector Type	Fascia Mounting Socket
Connector Size	Standard
Colour	Violet
Max. Temperature	120°C
+Positive Contact	Nickel Chromium
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS